

APAS[®] PharmaQC

Pilot secondary validation study results

July-23



Executive Summary – Pilot Secondary Validation study

Total 1,515 environmental monitoring plates collected during routine production

- 1,477 Environmental monitoring plates collected from Class A (including settle and finger dab plates)
- 38 contrived plates (11 pos, 27 neg) collected to increase positivity rate
- Data collected from LIMS system to provide reference result for each plate

Plates subsequently run through APAS Independence to assess performance using “real-world” data and identify any discrepant results between human and AI

- Measured growth vs no growth performance



Summary of performance

APAS PharmaQC demonstrates 100% sensitivity for microbial growth detection

Discrepant analysis investigations:

- APAS identified potential additional positive plate(s) missed by normal reading
- Top causes of disagreement between reference count and APAS:
 - Tape applied to plates during processing
 - Labels incorrectly applied to plates
 - Excessive text markings on the plate
 - Product identified on the plate




Results:

		APAS		
AZ		No growth	Growth	Total
	No growth	1390	110	1500
	Growth	0	15	15
	Total	1264	251	1515


- Positive percent agreement (PPA): 100% (15/15 correct)
- False negative rate (FNR): 0% (0/15 incorrect)
- Negative percent agreement (NPA): 92.6% (1390/1500 correct)
- False positive rate (FPR): 7.3% (110/1500 incorrect)

Example Results – Plate classifier result, single colony detection



APAS correctly identifies and counts single small colony on plate.

No interference from dish markings or labels



Example Results - APAS PharmaQC detects colonies missed by human reading

APAS INDEPENDENCE @Replay Mode Pathology Inc. About Sessions > Listing > Sample 1617/1877 → QC Sample ID Q

Signed in as operator, Sign out

▼ Report

Sample
010405612357632417230710102357041-
21002415544-910241

TSA 90mm, w/ Neut. +

Designation
Probable

Growth
1 colonies found.
wc: 1

Review reason
N/A

Count
1 colonies

▼ Image

Top-Lit Image

≡ # 🔍 🔍 ↻ 1:1 ⌛

shift+drag to measure clear


► Colonies

► Disc Panels

► Plate Details


► LIS Details

► About



Small white colony detected at edge of plate, missed during normal reading.

Example Results – Plate Counting Performance



APAS Count identifies 16 colonies which are confirmed by manual counting

Plate counting tools available in APAS Web-UI assist for human verification of results

Next Steps

Pilot primary validation study underway, testing performance and preliminary methods for formal testing. Testing scope aligned with Pharmacopeial requirements for:

- Linearity of Count
- Specificity
- Robustness
- Ruggedness
- Limit of Detection
- Precision
- Accuracy
- Operational Range
- Limit of Quantification
- Repeatability





For further information contact:

Jack Brown
Corporate Development Director
jack.brown@lbtinnovations.com

